





Cirrus

Model Number	PR-508M
Size	Medium
Weight Capacity	375 LB
	Maxi
Recline Positions	Comfort
Standard Fabrics ¹	6
Fabric Options	All
Chaise Seat	Yes
Coil Springs	Yes
Back Type	Biscuit
Head Rest Cover	Yes
Arm Covers	No
Battery Back Up	Yes
Overall Width	32"
Overall Height	43"
Width Between Arms	20"
Floor to Top of Seat	20"
Seat Depth	20.5"
Seat to Top of Back	29"
Seat to Top of Arm	5.5"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	23"
Floor to Top of Back Reclined	22"
Back Angle Closed	100°
Back Angle Reclined	170°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	26"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Removable Back Option Available	No
Smart Tek System Standard	Yes
Customizable *	No
Left Hand Control Option Available	Yes
Extra Pocket Option Available	Yes
Heat & Massage Option Available	Yes
Gel Memory Foam Option Available	Yes
Foot Rest Extension Option Available	Yes
Straight Lift Option Available	No
User Height (suggested)	5'4"-5'10"

¹ Buckskin Brisa is a standard, stock fabric on PR508M. Upcharge applies on all sizes.

* See Customize Chart for Dimension limits Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.

MAXICOMFORTTM CLOUD SERIES



Our patented technologies, combined with our superior kiln-dried northern maple hardwood frame design, provide the most comfortable lift and recline experience.



The Cirrus[™] features a combination of pocket coil seat springs surrounded by foam that molds to your lower body for the ultimate support and comfort.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.






Smart Tek

Our luxurious lift and recline chairs are equipped with electronics to help assure a secure lift every time you need it.



Literature is current at time of printing Golden Technologies reserves the right to make changes to the product or literature at any time.







Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com