



RESTORE
PEACE-OF-MIND
AND ACCESS
IN YOUR HOME

www.bruno.com





COMBINE QUALITY, COMFORT AND STYLE

When you want the highest quality stairlift combined with a stylized appearance, you want the Bruno Elite. The Elite offers worry-free access to all levels of your home. And the Elite's stylized design blends with your home's environment.

Let Bruno's Elite safely connect you to all levels of your home again.

- Made-in-the-USA quality
- Easy, reliable operation
- Smooth, quiet ride
- Luxury styling and comfort



Learn more at www.bruno.com


Even a Legend Can Use a Lift



“ We completed a sports room in our home to enjoy the football season with our family just before Bart had his strokes. Afterwards, it was not possible to use this area because Bart was unable to negotiate the stairs. A Bruno stairlift allows us to enjoy our entire home again. ”

Bart & Cherry Starr

“Green Bay Packers Hall of Fame Quarterback”



ELITE FEATURES

- Plush, generous-size seat with multiple height adjustments
- Adjustable width armrests
- Anodized rail with covered gear rack gives a sleek appearance
- Flip up arms, seat, and footrest create plenty of extra space on steps
- Adjustable footrest for maximum comfort
- Vertical rail design enables tight installation to within 6 in (152 mm) of wall
- Onboard audio/visual diagnostics for easy service
- Continuous charge strips along rail recharge unit wherever it stops
- Two wireless call/sends allow unit to be controlled remotely
- Two 12v batteries ensure dependable performance even in a power outage
- Soft start/stop gives a smooth ride from start to finish




Ask your local dealer about upholstery and color options.



Offset swivel seat makes it easy to get on/off.

“The Bruno stairlift has eliminated an enormous amount of anxiety and stress. We feel grateful for it every day.”

A. Wambold - VA




Close installation to wall.

ELITE OPTIONS



Power folding footrest automatically flips up/down when seat is raised/ lowered.



Power or manual folding rails for narrow hallway or when doorway is at the bottom of stairs. Manual operation or push button automatic.



Power swivel seat for effortless exit. Controlled on chair arm or wireless call/send.



Larger seat and footrest for increased comfort.



"The stair lift is great, very solid and looks good too."

L. UHRING - MI



TECHNICAL SPECIFICATIONS

Standard Seat Size: 17.5 in wide x 15 in deep
(445 mm wide x 381 mm deep)

Larger Seat Option: 21.5 in wide x 16 in deep
(546 mm wide x 406 mm deep)

Standard Footrest Size: 10 in long x 9 in wide
(254 mm long x 229 mm wide)

Larger Footrest Option: 12 in long x 11 in wide
(305 mm long x 279 mm wide)

About Bruno

Bruno Independent Living Aids has helped improve the lives of people with limited mobility for over 30 years. Engineered and manufactured in the USA, Bruno stairlifts and vertical platform lifts help people better access homes, public buildings and businesses, and Bruno's full line of scooter/powerchair lifts help people easily transport their mobility devices. In addition, Bruno's Valet Signature Seating is the world leader in revolutionary turning automotive seats to assist people in getting in and out of their vehicles more easily. A family-owned business, Bruno takes pride in their unwavering focus on quality and safety and their position as industry innovation leader.

Proudly made in the USA.

 **LINCOLN
MOBILITY**

402-421-8800 | www.LincolnMobility.com



 **BRUNO®**

Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066

www.bruno.com

© Bruno Independent Living Aids, Inc. 2018
2018-01-16-Elite-6-PageV1.indd • MKT-SL-2