



Bruno Residential
Vertical Platform Lift

TAKE YOURSELF
HIGHER





RESIDENTIAL **VPL**

FIT YOUR SPACE AND LIFE

MINIMUM SPACE. MAXIMUM REACH

Space-efficient Bruno wheelchair lifts can access up to 14 feet.

INDOOR OR OUTDOOR

Weather-resistant controls and specialized coating protect the unit. Combat extreme winters with the optional cold weather package.

HEAVY-DUTY LIFT CAPACITY

Lift capacity of 750 lbs (136 kg).

2-YEAR LIMITED WARRANTY

www.bruno.com



Give Your Life a Lift

Gain access to your porch, deck, or floor-to-floor without sacrificing space with a compact Bruno residential vertical platform lift. The quiet, smooth lift accommodates any wheelchair, power chair or scooter. Give yourself the freedom to get in and out of your home with just a push of a button. You're in control.

COUNT ON BRUNO QUALITY

EASY-TO-USE – Pressure control rocker switch up/down operation. Automatic self-lowering ramp folds down for smooth roll on/off access; folds up when in use to create a safety barrier.

HIGH-PERFORMANCE – Lift capacity of 750 lb/340 kg. Smooth, quiet ride. Made-in-the-USA quality.

SAFE – Unit shuts down when obstruction is touched under platform. Emergency stop switch on platform.

RELIABLE – Continuously charged DC battery-powered motor on lifts 10 ft and above ensures operation even in a power outage. Affordable AC in 4- and 6 ft units (DC battery powered optional)
Gold warranty: 2 years, major components; 1 year, parts.



OPTIONS TO MEET YOUR NEEDS

- Top landing gate
- Landing call/send station
- Platform gate
- Multiple platform sizes and configurations
- Cold-weather package
- Paddle controls
- DC battery powered in 4- and 6-foot units





About Bruno

Bruno Independent Living Aids has helped improve the lives of people with limited mobility for nearly over 35 years. Engineered and manufactured in the USA, Bruno stairlifts and vertical platform lifts help people better access homes, public buildings and businesses, and Bruno's full line of scooter/powerchair lifts help people easily transport their mobility devices. A family-owned business, Bruno takes pride in their unwavering focus on quality and safety and their position as industry innovation leader.

Proudly made in the USA.




LINCOLN MOBILITY

402-421-8800 | www.LincolnMobility.com



CONSUMERAFFAIRS



BRUNO®

Bruno Independent Living Aids, Inc.
1780 Executive Drive, Oconomowoc, WI 53066

www.bruno.com

© Bruno Independent Living Aids, Inc. 2019
2019-03-07-RES-VPL-4pg.indd • MKT-VPL-TT